



## Colour (Color) Standards

ASTM, Saybolt, Platinum-Cobalt (Hazen), Gardner,  
European and United States Pharmacopoeia

Reagecon manufactures the full range of ASTM, Saybolt, Platinum-Cobalt (Hazen), Gardner, European Pharmacopoeia and United States Pharmacopoeia Colour Standards for use with ASTM, APHA, ACS, European and United States Pharmacopoeia standard methods. The ASTM standard methods include D1500, D6045 and D1209. The products can be used to calibrate, control, qualify and validate colour measurement instruments.

- ASTM Colour Standard Sample A05- A7
  - Saybolt Colour Standards S +30 to S -15
  - Platinum-Cobalt Scale No. 0 - No. 1000
  - Gardner Colour Standards GARD01 - GARD18
  - European Pharmacopoeia Standards (Opalescence, Primary and Standard Solutions)
  - United States Pharmacopoeia Standards
  - Chinese Pharmacopoeia Standards
- 
- ✓ Produced in accordance with ASTM (D1500, D6045, D1209) APHA, ACS, EP & USP methods
  - ✓ Consistency of product - Independent, Traceable, Certified
  - ✓ For use as Calibration and/or Quality Control Standards
  - ✓ Presented in high quality tamper evident bottles
  - ✓ Various pack sizes available
  - ✓ Ready to Use
  - ✓ Customised Standards available
  - ✓ Certificates of Analysis and Safety Data Sheets available online

