

Redox Standards

Reagecon produces a complete range of high-quality Redox Standards for the control, validation or qualification of Redox measurement. In fact, our range are the widest available on the market, have a very high specification of ± 5 mV and guaranteed stability throughout their entire shelf life, even after opening the bottle.

Products come with comprehensive documentation that includes Safety Data Sheets and detailed Certificates of Analysis with proven verifiable accuracy and fully budgeted Uncertainty of Measurement.

The products are available in several different pack sizes that include 500ml, 10 litre and 10 litre bag-in-box options.

Customised products maybe available upon request.



- Extensive range 124 mV – 650 mV @ 25°C
- Very high specifications (± 5 mV)
- Extensive technical advice on the measurement techniques available
- All products certified with proven verifiable accuracy and uncertainty of measurement
- Widest range of values and pack options available in the market
- Ready to Use
- Certificates of Analysis and Safety Data Sheets available online

For further information or enquiries, please contact:

Tel: +353 61 472622 | Email: sales@reagecon.ie

www.reagecon.com