

Brix / Refractive Index Standards

Reagecon manufactures several ranges of Brix / Refractive Index (RI) Standards for ease of use when controlling all types of Refractometers, irrespective of brand.

- ✓ ICUMSA compliant Brix Standards (Sucrose)
- ✓ Non ICUMSA compliant Stabilised Brix Standards (Sucrose)
- ✓ RI Standards manufactured from Sucrose
- ✓ RI Standards manufactured from Solvents
- ✓ Sucrose in Water (Brix) ISO 17034 Certified Reference Materials (see separate flyer)

Stabilised Brix Standards have an extended shelf life of 1 year and Refractive Index Solvent based Standards offer a 2 year shelf life.

Products represent excellent value for users that are not required to follow ICUMSA Guidelines.

These product ranges are used primarily either as a Calibrant or Analytical Quality Control Solution.



- ✓ Extended Shelf Life of:
 - 1 year Stabilised Sucrose products
 - 2 years Solvent Based products
- ✓ Extensive range 0 – 67.5% Brix
1.33310 – 1.65812 nD Refractive Index
- ✓ High accuracy:
 - ± 0.11% @ 20°C Brix
 - ± 0.00014 nD units Refractive Index Stabilised Sucrose
- ✓ ISO 17025 Accredited for 5 - 60% Brix
and 1.3310 – 165812 nD Refractive Index
- ✓ Consistency of product-Independent, traceable, certified
- ✓ Presented in a convenient high quality dropper bottle
- ✓ Can be used with any brand of Refractometer
- ✓ Ready to Use
- ✓ Certificates of Analysis and SDS Sheets available online

For further information or enquiries, please contact:

Tel: +353 61 472622 | Email: sales@reagecon.ie

www.reagecon.com