

A CALIBRE SCIENTIFIC COMPANY

Sucrose in Water (Brix) ISO 17034

Reagecon manufacture a range of Sucrose in Water (Brix) Certified Reference Materials (CRMs), which are manufactured and certified in accordance with the requirements of ISO 17034. These products are used primarily either as a Calibrant or Analytical Control Solution in Refractive Index based methods of Brix value determinations, they can also be used to validate appropriate test methods or qualify a refractometer for use in a regulated industry. Products have a shelf life of 20 weeks and are produced in accordance with ICUMSA guidelines.



- Extended Shelf Life (20 weeks)
- In accordance with ICUMSA guidelines
- Extensive range 1 60% mass/mass
 Sucrose in Water solutions
- Product measurement uncertainty is computated on a batch to batch basis, guaranteed to never exceed ± 0.15 °Brix
- Consistency of product Independent, Traceable, Certified

- Presented in a convenient high quality dropper bottle
- Can be used with any brand of Refractometer
- Ready to Use
- Certificates of Analysis and SDS sheets available online

For further information or enquiries, please contact: Tel: +353 61 472622 | Email: sales@reagecon.ie www.reagecon.com