

Antimony pH Buffer Solutions

Reagecon is a world class manufacturer of pH Buffers, offering the most comprehensive range available to the marketplace, designed to suit all requirements

In addition to conventional pH Buffer Standards, Reagecon manufacture a range of specialist Buffer Standards formulated for use in intra-gastric pH analysis using Antimony or equivalent electrodes

- ✓ pH values 1.07, 4.00 & 7.01 at 25°C
- ✓ High Specification ± 0.05 pH units
- ✓ Phosphate Free
- ✓ For use with pH testing systems and Antimony or equivalent probes
- ✓ Consistency of product, Independent, Traceable, Certified



QUALITY PRODUCTION

pH Testing is ISO17025 Accredited
Values specified at 25°C

Manufactured to exacting specifications with an extended shelf life & stability

TRACEABILITY

Directly traceable to the IUPAC pH scale by an unbroken chain of traceability. Reagecon achieve this traceability through a series of comparisons, with the key reference materials being Standard Reference Materials (SRMs) manufactured by NIST

STABILITY

Guaranteed stability throughout entire shelf life, even after opening the bottle, eliminates the requirement to open a fresh bottle of standard every time it is used

LABELLING

Temperature dependence data is printed on the label as are lot numbers and expiry dates for user convenience

COA & SDS

Certificates of Analysis and Safety Data Sheets available online

For further information or enquiries, please contact:

Tel: +353 61 472622 | Email: sales@reagecon.ie

www.reagecon.com