

## *pH Buffers ISO 17034 and ISO 17025*

Reagecon manufactures the most comprehensive range of pH Buffer Solutions available to the marketplace, which have been designed to suit all requirements.

- ✓ Laboratory Grade Buffers, pH range 1-13
- ✓ Colour Coded Buffers for Calibration pH 4, 7, 10
- ✓ Professional Range (Technical, NIST/DIN & High Resolution Buffers)
- ✓ Low Ionic Strength Buffers
- ✓ Borax Free
- ✓ pH Buffer Capsules
- ✓ Various packaging options for ease of use e.g. Colour Coded, Twin Neck Bottle and Bag-in Box

### Electrode Care & Maintenance

Reagecon manufactures a complete range of products designed to clean, refill, regenerate and store your electrodes.



### QUALITY PRODUCTION

Manufactured and certified in accordance with the requirements of ISO 17034 or ISO 17025 at specific pH values and temperatures (e.g. 20°C/25°C). Both our ISO 17034 pH Buffer Certified Reference Materials (CRMs) and ISO 17025 pH Buffer Solutions are manufactured to exacting specifications with extended shelf life & stability

### TRACEABILITY

Directly traceable to the IUPAC pH scale by an unbroken chain of traceability. Reagecon achieve this traceability through a series of comparisons, with the key reference materials being Standard Reference Materials (SRMs) manufactured by NIST

### LABELLING

Temperature dependence data is available via the product label as are lot numbers and expiry dates for user convenience

### COA & SDS

Certificates of Analysis and Safety Data Sheets available online

For further information or enquiries, please contact:

Tel: +353 61 472622 | Email: [sales@reagecon.ie](mailto:sales@reagecon.ie)

[www.reagecon.com](http://www.reagecon.com)