



Reagecon

A CALIBRE SCIENTIFIC COMPANY

November 2024

Reagecon Cryoscope Standards

Dear Customer,

Reagecon's Cryoscope Standards are traceable to the National Institute of Standards and Technology (NIST) and are manufactured and certified for use on all Cryoscopes compliant to International Reference Standard ISO5764 / IDF108 "Milk – Determination of freezing point – Thermistor cryoscope method."

If you have any queries regarding Reagecon's Cryoscope Standards, please do not hesitate to contact us at quality@reagecon.ie

We thank you for your continued business.

Regards

Darren McGrath
Quality Director

Quality Department
Reagecon